

# Competent Cell Preparation and Transformation

## [Preparing the competent cells]

### Reagent

1. TSS (Transformation and Storage Solution for chemical transformation) /50 mL
  - 85% LB medium : 1.25 g x 0.85 = 1.0625 g
  - 10%(wt/vol) Polyethylene glycol (MW 8000) : 5 g
  - 5% DMSO (vol/vol) : 2.5 mL
  - 50mM MgCl<sub>2</sub> anhydrous (pH 6.5) : 0.238 g

Autoclave in a **glass bottle** or filter sterilize. Store at 4 °C for < 2 weeks

2. **2x LB** medium /200 mL

- LB medium : 6.2 g
- NaCl<sub>2</sub> anhydrous : 3.8 g

Autoclave in a 500-mL **baffled flask**.

### Methods

1. Streak on a LB plate (added antibiotic if cells have antibiotic resistant). Incubate the plate at 37 °C overnight.
2. Pick a single, well-isolated colony and inoculate it into 5 mL of LB broth (plus antibiotic). Incubate at 37 °C overnight with shaking at 200 rpm.
3. Transfer 2 mL of the saturated overnight culture to a sterile 500-mL flask containing 200 mL of 2x LB medium (**do not add antibiotic although cells have antibiotic resistance**). Incubate the cells at 37 °C with the shaking at 200 rpm, until OD<sub>600</sub> reach **0.6 ~ 0.8**.
4. When the culture reaches an OD<sub>600</sub> of 0.8, chill the flask on the ice for 20 min (completely cool down the culture) and then collect the cells by centrifugation at 1,500 rpm for 5 min at 4 °C.
5. Resuspend the cells in 5 mL of ice-cold TSS solution. Now the competent cells are ready to be transformed.
6. Aliquot 110 uL competent cells to ice-colded 1.5 mL tube. If they are not immediately used, cells can be stored at 4 °C for maximum of 6 hours without significant loss of competency. The same competent cells can also be stored at -70 °C for long-term storage (pre-treat with liquid nitrogen).

## [Transforming the cells]

1. Thaw -70 °C competent cells first on ice. Add DNA (**do not exceed 10% vol. of competent cells**) to ice cold 50 uL competent cells.
2. Incubate on the ice for 20 ~ 30 min.
3. Heat shock at 42 °C for 1 min.
4. After heat shock, put on ice for 2 min.
5. Add 200 uL of **pre-warmed** LB broth.
6. Shake and incubate at 37 °C for 60 min.
7. Plate on the **pre-warmed** appropriate agar plates and incubate plates at 37 °C, overnight.